

## ***DAPI staining of *S. pombe* cells***

### ***Fixation***

#### ***Glutaraldehyde fixation (beautiful way)***

***Nucleus, cell wall, and septum can be detected***

Culture cells in YES until OD600=0.5~1.0, @30°C or appropriate temperature, 25ml

Get 900 ul of culture in a microfuge tube.

Add 100ul of glutaraldehyde.

Immediately vortex, for 30sec.

Place the tube on ice for 30min.

Wash cells in ice-cold 1xPBS three times

(Centrifuge, discard sup, add 1ml of 1xPBS to the pellet, vortex)

Suspend cells in 100 ul of 1 x PBS, keep tubes on ice

**\*\* Cells can be stored @4°C (take pictures within 24hr, but the sooner the better)**

Centrifuge briefly to collect cells at the bottom of tubes, keep tubes on ice

Mix 1ul of cells (from the bottom) and 1ul of DAPI (50ug/ml) on a slide glass.

**\*\* in the microscope room.**

Place a cover slip (22x22) over the cells

#### ***Ethanol Fixation (dirty way)***

***Nucleus can be detected***

10ml of cell culture in 15ml conical tube

Centrifuge, 2400rpm, 4°C, 3min

Discard sup

Add 70% Ethanol to the cell pellet

Vortex well, incubate cells @RT for 10min

**\*\*Cells can be stored @-20°C**

Wash cells in ice-cold 1xPBS three times

(Centrifuge, discard sup, add 1ml of 1xPBS to the pellet, vortex)

Suspend cells in 100ul of 1xPBS, keep tubes on ice

**\*\* Cells can be stored @4°C (take pictures within 24hr, but the sooner the better)**

Centrifuge briefly to collect cells at the bottom of tubes, keep tubes on ice

Mix 1ul of cells (from the bottom) and 1ul of DAPI (50ug/ml) on a slide glass.

**\*\* in the microscope room.**

Place a cover slip (22x22) over the cells

#### ***No Fixation (the most beautiful way)***

***Nucleus, cell wall, and septum can be detected***

1.5ml of cell culture in a microfuge tube

Wash cells in ice-cold 1xPBS **once**

(Centrifuge, discard sup, add 1ml of 1xPBS to the pellet, vortex)

Suspend cells in 100ul of 1xPBS

Centrifuge briefly to collect cells at the bottom of tubes, keep tubes on ice

\*\*\* **Cells cannot be stored. Take all pictures within 30min.**

Mix 1ul of cells (from the bottom) and 1ul of DAPI (50ug/ml) on a slide glass.

\*\* in the microscope room.

Place a cover slip (22x22) over the cells