


E. coli Transformation

by Chiaki Noguchi



If your plasmid has an ampicillin resistant gene, spread your E. coli cells on to LB + 2 mg/ml ampicillin plate.

We have 50 mg/ml ampicillin in 50% EtOH at -20°C .

For one LB+amp plate:

Take $50 \mu\text{l}$ of this stock solution into $50 \mu\text{l}$ of sterilized milliQ water in a 1.5 ml tube to make a dilution mixture.

(If you have more plates, make more of this mixture.)

And spread $100 \mu\text{l}$ of it to a LB plate with 6~8 glass beads in.

Making Competent Cells

by Chiaki Noguchi

day 1

17:00 Spread DH5a on LB plate from glycerol stock

Place the plate @37°C over night

day 2

15:00 single colony 1 }
LB 5ml } x 2 sets in 50ml tubes

Shake @33°C over night

>>>Place 2 sets of 100ml LB in 500ml flasks @33°C room


day 3

10:00 Put 1 tube of culture cells into warmed 100ml LB in a flask


Shake @33°C 3.5hr

>>>From here, keep it cold!!<<<

13:30 Put 200ml culture cells into 50ml tubes (x4)
on ice 5min

 3500rpm 5min @4°C

Suspend cells with 20ml Tfb1 each tubes

 3500rpm 3min @4°C

Suspend cells of 4 tubes with 8ml Tfb2

Aliquot 100ml in 0.6ml tubes

Put them into liquid nitrogen

Keep it @-80°C

Transformation Buffer 1 & 2 for competent cells

by Chiaki Noguchi

>>> Tfb1

	stock			
30mM KoAC	4M	1.88	ml	(4M...19.63g/50ml)
0.1M RbCl		3.02	g	
10mM CaCl ₂	1M	2.50	ml	
50mM MnCl ₂		2.47	g	
15% Glycerol	60%	62.50	ml	
DW		to 250.00	ml	

↓
0.22mm filtration (40ml/tube x6)

↓
@ -20°C

>>> Tbf2

	stock			
10mM MOPS		0.52	g	
10mM RbCl		0.30	g	
75mM CaCl ₂		2.08	g	
15% Glycerol	60%	62.50	ml	
DW		to 250.00	ml	

↓
pH~4.7 to 6.5 (by KOH)

↓
0.22mm filtration (8ml/tube x6 & 50ml/tube x4)

↓
@ -20°C