

## **Immunofluorescence Microscopy of *S. pombe* cells.**

### **Day1:**

Cell culture in YES (10ml), 30°C. (25°C for ts cell lines)

### **Day2:**

If cells are overgrown (OD: more than 1.5 in YES, 1.0 in EMM), dilute cells in medium at OD=0.2 and culture cells for two generation times (4h in YES, 6h in EMM; if cells grow like wild-type). Optimal condition for microscopic analysis is an OD of 0.5~1.0 in YES, 0.4~0.8 in EMM.

#### **Fixation**

Add 1ml of formaldehyde (30%) to the cell culture (10ml)

Vortex well and wait for 45sec @ RT

Add 100µl of glutaraldehyde (25%)

Vortex well and incubate for 1h @ 30°C (@25°C for ts) in a shaker

Centrifuge, 2000rpm, 3min, remove supernatant completely

Suspend cells in 1ml of PEM and transfer the suspension to 1.5ml tube

Wash cells in 1ml of PEM for three times

{Centrifuge (10krpm, 30sec), remove sup, add 1ml of PEM}

At this point you have cells in 1ml of PEMS

#### **Cell wall digestion, permeabilization**

Add 5ul of 200 x Zymolyase, vortex briefly to mix

Incubate, 70min, @37°C, invert tube occasionally

Check cells if cell wall is digested under the microscope

Centrifuge, 8krpm, 10sec

Add 200µl of PEMST, suspend cells gently, wait for 2min

Wash cells in PEM for three times

{centrifuge (12krpm, 60sec), remove sup, add 200µl of PEM}

#### **Quenching**

Centrifuge, 12krpm, 60sec, remove sup

Suspend cells in 200µl of PEM containing 1mg/ml Sodium borohydride

Incubate, 5min, RT, with the lid open

Repeat the quenching steps another two times

Preparation of Sodium borohydride solution

Stand by 1mg of Sodium borohydride in a 1.5ml tube (prepare three tubes)

Add 1ml of PEM just before quenching

Use 200µl to suspend cells for one quenching reaction

Wash cells in 200µl of PEM for three times

{centrifuge (12krpm, 60sec), remove sup, add 200µl of PEM}

#### **Blocking, antibody reaction**

Centrifuge (12krpm, 60sec), remove sup, suspend cells in 200µl of PEMBALG

Centrifuge (12krpm, 60sec), remove sup, suspend cells in 100µl of PEMBALG  
Incubate for 1h @RT on the nutator  
Centrifuge (12krpm, 60sec), suspend cells in 100µl of PEMBALG with the 1<sup>st</sup> antibody  
Example: Anti-Myc (9E10), 1/2000  
Incubate o/n @RT on the nutator

### **Day3:**

Wash cells in 100µl of PEMBALG for three times  
{centrifuge (12krpm, 60sec), remove sup, add 100µl of PEMBALG}  
Centrifuge (12krpm, 60sec), remove sup, suspend cells in 100µl with 2<sup>nd</sup> antibody  
Example: 1/200  
Cover the tube with aluminum foil (2<sup>nd</sup> antibody is light sensitive)  
Incubate o/n, RT, on the nutator

### **Day4:**

Wash cells in 100µl of PEMBALG for three times  
{centrifuge (12krpm, 60sec), remove sup, add 100µl of PEMBALG}  
Centrifuge (12krpm, 60sec), remove sup, suspend cells in 50µ of PBS+0.2µg/ml DAPI  
Immediately, Centrifuge (12krpm, 60sec), remove sup  
Suspend cells in 20µl of PBS+0.02% Azide, store cells at 4°C  
You are now ready to go to the microscope room