

Kunkel Method

By Chiaki Noguchi

Modified by Jordan Rapp (1/20/07)

= ssDNA prep =







day 1

- 1) Transform RZ1032 with plasmid DNA on (amp) plate @37°C

day 2

- 2) Culture 1ml LB(amp) with a single colony @37 °C for about 5.5hr
- 3) Take 0.2ml RZ-cell to 1.5ml tube, mix 0.1ml helperphage (promega/R408)
Incubate for 30 min at 37 °C
- 4) Transfer all to 20ml LB in 50ml flask, culture O/N @37 °C
w/ 40µl 50mg/ml amp + 10µl 10mg/ml uridine

day 3

- 5) Transfer culture to polypropylene bottle
 12000rpm 15min RT
- 6) Transfer sup to new bottle
 12000rpm 15min RT
- 7) Transfer sup to new bottle again (glass tube)
add 5ml 2.5M NaCl / 20% PEG8000, on ice 30min
 9000rpm 20min 4 °C
- 8) Discard sup
 9000rpm 20min 4 °C
- 9) Discard sup completely using pipetman
- 10) Suspend with 1ml TE, transfer to two 1.5ml tubes
phe x 2times, chl x 1time
- 11) Add 1/10vol. 3M NaoAc + EtOH @ -80 °C 30min
 Max 10min 4 °C
- 12) Suspend with 0.2ml TE, merge two to one 1.5ml tube
phe/chl x 1time, chl x 1time
EtOH, 70% EtOH wash

- 13) Suspend with 25µl TE
- 14) Check 1µl on TAE agarose gel, keep -20 °C
(good for quite long time)

by Chiaki Noguchi

= Mutagenesis =

Kination

5pmol/ μ l oligo	15 μ l
10x Kination buffer	2 μ l
10mM ATP	1 μ l
<u>T4-polynucleotide kinase</u>	<u>2 μl</u>
	20 μ l

37oC 30min

65oC 20min ==> P-oligo (Keep @ -20oC)

Annealing

UssDNA	1 μ l	1 μ l
P-oligo	1 μ l	-
10x Annealing buffer	1 μ l	1 μ l
<u>D.W.</u>	<u>7 μl</u>	<u>8 μl</u>
	10 μ l	10 μ l

70oC 10min

Cool down to RT for more than 3 hours, then put it on ice
(DON'T RUSH! This is important!!)

Synthesis

add 10x Synthesis buffer 1 μ l
T4 DNA ligase (3units/reaction)
T4 DNA polymerase (1unit/reaction)

0 °C 5min

25 °C 5min

37 °C 90min

add STOP buffer 90 μ l

Chcek 10 μ l by agarose gel electrophoresis

Transformation E.coli with 10 μ l of product

= Reagents =

10x Kination buffer (2 μ l/reaction)

700mM Tris-HCl pH7.4

100mM $MgCl_2$

50mM DTT

10x Annealing Buffer (1 μ l/reaction)

200mM Tris-HCl pH7.4

20mM $MgCl_2$

500mM NaCl

10x Synthesis buffer (1 μ l/reaction)

5mM each dNTP

10mM ATP

100mM Tris-HCl pH7.4

50mM $MgCl_2$

20mM DTT

T4 DNA ligase (3 units/reaction)

T4 DNA polymerase (1 unit/reaction)

STOP buffer (90 μ l/reaction)

10mM Tris-HCl pH7.4

15mM EDTA

by Chiaki Noguchi