

SDS-PAGE

(Poly Acrylamide Gel Electrophoresis)

Reagents:

Laemmli Buffer

- 1) Clean a set of glass plates (the one is flat, the other has spacer) with 100% Ethanol. Use Kimwipe to clean plates. Dry them up for few seconds.
- 2) Set glass plates as described in the Bio-Rad direction. The side you clean is inside.
- 3) Prepare the Lower (Separating) Gel solution (See recipe). Pour the gel solution between the glass plates. Cover the gel surface with 0.25ml of 2-propanol. Wait for 20 to 30 min until the gel becomes solidified.
- 4) Discard the 2-propanol layer. Wash the gel surface with MilliQ water, twice. Remove all the remaining water with a piece of 3MM paper. But don't touch the surface of Lower Gel.
- 5) Prepare Upper (Stacking) Gel solution (see recipe). Pour with the upper gel solution on top of the solidified Lower gel with a comb in between the plates. Wait for 20 to 30 min until the gel becomes solidified.
- 6) Set the gel sandwich into Bio-Rad Mini-PROTEAN 3 Cell system. Fill the chamber with Laemmli Buffer. Don't mix inner buffer and outer buffer.
- 7) Prepare your protein samples. Mix well and centrifuge to collect the liquid at the bottom of the tube. Thaw the Pre-Stained Protein Ladder Marker. Mix well and centrifuge to collect the liquid at the bottom of the tube.
- 8) Apply marker and samples with Hamilton Syringe (15ul for a large gel, 8ul for a mini gel). The applying volume of your samples are depend on a case.
- 9) Put the lid on. Start running the gel.
 - 35 mA constant current for 1 gel for about 30 min
 - 50 mA constant current for 2 gel (one container) for about 30 to 40 min
 - or
 - You can reduce current to run the gel slowly. For example:
 - 25 mA constant current for 1 hour
- 10) Run the gel until the front line reaches at the bottom.