

Pulsed-Field Gel Electrophoresis of Mammalian Genomic for DNA Damage Determination

Reagents:

Standard Tissue Culture conditions, equipment, and protocols.

BioRad Low Melt Agarose

BioRad MegaBase Agarose

BioRad PFGE CHEF-DRII System

BioRad Plug Molds

Digestion Buffer:

1% Sodium Laurel Sarcosine

0.2% Sodium Deoxycholate

1 mg/mL Proteinase K

100 mM EDTA, pH 8.0

1x Wash Buffer:

20 mM Tris-HCl, pH 8.0

50 mM EDTA, pH 8.0

Protocol:

- (1) Harvest Cells via Trypsinization, first rinse cells briefly in 1xPBS, trypsinize and quench with media. Pellet cells, resuspend in 1x PBS.
- (2) Count cells with a hemacytometer.
- (3) Pellet cells again, flash freeze pellet in Dry Ice/EtOH, transfer to -80°C freezer.
- (4) Resuspend pellet in 2% Low Melt Agarose in 0.5xTBE buffer, resuspend each sample so as to equalize cell concentration in the LMT solution, mix well.
- (5) While still liquid transfer Cell/LMT solution to Mold Plugs, taking care not to introduce bubbles to plugs.
- (6) Solidify plugs by placing on ice.
- (7) Lyse cells in plugs, by incubating in Digestion buffer at 55°C for more than 24hrs.
- (8) Wash plugs in Wash Buffer, 4x.
- (9) Store plugs in 0.1x Wash Buffer.
- (10) Prepare 0.7% gel in BioRad MegaBase agarose in 0.5x TBE. Insert plugs, and seal with LMT, so they are stable inside wells. NO BUBBLES.
- (11) Add a sufficient volume of 0.5x TBE to the CHEF-DRII apparatus and set to 14°C. Pump strength of 70%.
- (12) Transfer gel to CHEF-DRII apparatus.
- (13) Run first stage of Electrophoresis: 30 hours, pulse gradient from 30s to 120s, and a field strength of 1.9V/cm³.
- (14) Run second stage of Electrophoresis: 51 hours, pulse gradient from 120s to 42 minutes, and a field strength of 1.9V/cm³.
- (15) Remove gel from apparatus.
- (16) Stain gel in EtBr (0.5 µg/mL) in dH₂O for 30 minutes.
- (17) Destain gel in dH₂O for 3 hours.
- (18) Image gel